


Case Submission Using VetPacs

Click on the “Search” button under the Patient search heading. A timeframe can be selected by using the drop down arrow to help narrow the search. Double click on the patient study from the “Archive” to select the patient study. If the history has not already been entered, add the history and click the “Save Exam” button.



If this is a new patient, click the “Create/Edit Exam “ button. Fill in all of the patient information and history and click the “Save Exam” button.

Create Exam

Exams:
Clinic Name: DarkHorse Demonstration Center Exam Date: 10/31/2006 12:31

Patient Info:
Name: BRANDON Species: CANINE Breed: MIXED
Sex: MALE-NEUTERED Weight: Client First: Client Last: BRANDON Patient ID: 69019
Birth Date: Age: Yrs. Mos.

Exam Info:
Reference ID: Accession ID: Category:

Presenting Problem:
Notes to Specialist - Patient history, Laboratory results, Clinical signs (Does not appear on written report)


Veterinarian Info:
Veterinarians: Add Veterinarian Edt Veterinarian

Macros:
Export Macros Manage

Sent Exam Info:
Specialist: Invoice Number: Send Date:

Save Exam Acquire Images Cancel

Select the “Process to the Outbox” button. The type of study can be selected in the “Categories” heading by clicking on the appropriate choice. The “Radiographs Only” button will default to a 24-hour turnaround for the final report. The “Stat Radiographs” button will set the consult on a 2-hour turnaround. To change the consult to 6 or 12 hour turnaround, contact the DarkHorse office. The “Store” key submits the study to the PACS system with no report. Under the heading “Telemedicine Web Service URL’S” select DarkHorse as the destination for the study. Be sure to click the “Select All Images” button. Then click on the “Process Exam to Outbox” button.



From the Outbox select the exam and then click on the “Send Exam” button. To process multiple exams at the same time, send all of the exams to the outbox and select the “Mark All” button followed by the “Send Exam” button.

VetPACS - Outbox

File Configuration Exams Window Help

Archive Outbox Inbox Search

SEND EXAM EDIT EXAM MARK ALL REMOVE EXAM TO ARCHIVE

Exams							
Exam Date	Client Name	Patient Name	Veterinarian	Clinic	Category	Status	# to Send
11/11/2005	BRANDON,	BRANDON	User, Demo	DarkHorse Demonstration Center	RADIOGRAPHE ONLY	NEW/OK	7 of 7