HOW TO LINK STUDIES TO ACTIVE CONSULTS ON AIS

** If you have additional studies that you would like interpreted with your In Progress or Awaiting Media consult, here is how to do it:

1. Locate the study with a submitted consult you would like to link by clicking the ► button adjacent to the patient name.



This will expand the study so you can see all studies for the given patient.

2. Then use the "Select Action" drop down menu for the images without an associated consult to "LINK STUDY TO A CASE".



4. When the new screen pops up, you can simply check the circle and then hit the submit button

Lir	nk Study For		_						
Pat	ient Info From Selected Studie	es							
	Case Number: N/A Patient: TEST Owner: TEST, TEST Species: Canine Breed: Dachshund DOB: N/A Sex: U		(TEST	TEST	2012/08/09	SC 1	
	S	elect appro	priate	consult for linkin	g with selected	studies			
	Patient Info From Consultation			n	Study Type	Study [Date U	rgency	
	Case Number: N/A Patient: TEST Owner: TEST, TEST Species: Canine Breed: Dachshund DOB: N/A Sex: U Matches patient for sele	cted study		Ra	diology	Thu, 9 Aug	2012 3 hr	S	
Click Sub	mit.								
	Subm	it	ance	l and Retur	n To PACS				
						Deve			
NOW YOU Displaying 1 -	can see that both	studie August 09, 2012	s are 2 9:41 am	e linked:					
	Actions	Patient/Owner	Species	Hospital/Clinic	Submitted Date	Required Date	Time Remaining	Study Type	1
• • 🖳	() Select Action +	TEST/TEST	Canine	Antech Imaging Services Demonstration Center	2012/08/09 9:00 AM PDT	2012/08/09 12:00 PM PDT	2 hours, 18 minutes Remaining	Radiology	
		TEST/TEST	Canine	Services Demonstration Center	2012/08/09 9:00 AM PDT	2012/08/09 12:00 PM PDT	2 hours, 18 minutes Remaining	Radiology	

Powered by DarkHorse Medical Ventures Copyright © by Dark Horse Medical Ventures. All Rights Reserved